

DERMATUFF Sleeve Dimensions for the assistance of customers

This Information is useful particularly if you are considering wearing the sleeves on your legs rather than wearing socks with feet.

It should be noted that these Dimensions are minimum and not exact. The sleeves are made on a sock machine and because of the nature of the machine and of the materials there can be variation in the sizes produced.

There is a considerable degree of stretchiness in either direction in order to allow for the products to fit a greater variety of actual limb sizes.

Some slight pressure is necessary in order to the sleeves not to slip down but it is important not to apply to much pressure to the legs if you have vascular problems, the pressure applied increases with the width expansion but also the length decreases. If you have concerns about this then please refer to us on 01752 880324, or your medical advisor.

LARGE

Circumference at top: 42 cm

Circumference at bottom: 32 cm

Length top to bottom: 45 cm

MEDIUM

Circumference at top: 36 cm

Circumference at bottom: 24 cm

Length top to bottom: 37 cm

SMALL

Circumference at top: 30 cm

Circumference at bottom: 24 cm

Length top to bottom: 32 cm